

Shedding fresh light on particle characterisation

Malvern Instruments, Malvern, UK, has released a new version of software for its Morphologi G3 particle characterisation system. This fully automated image analysis-based technology delivers particle size and shape information for wet and dry samples, as well as for the analysis of foreign particulate matter collected on filters. With the new software, users can now more closely define their own light settings in the standard operating procedures, allowing for darkfield and polarizer measurements. The Morphologi G3 automated particle characterisation system provides statistically significant, high-quality image data for thousands of particles in one measurement. www.malvern.com/morphologi



Particle image from the Morphologi 3